



## MALDI Biotyper

- Fast & Accurate Identification of Microorganisms



MALDI Biotyper: an alternative  
to identify microorganisms

Dr Delphine Vincent  
11 Nov 2013

# Introduction

## The MALDI Biotyper System:

- Highly Accurate
- Applicable to a Wide Range of Microorganisms
- Much Faster than Traditional Methods
- Cost Effective
- Robust and Easy to Use

- Traditional methods of identifying microorganisms
  - Gram stain
  - Culture on selective media
  - Biochemical tests
  - 16S sequencing
- Some limitations of traditional methods
  - Time consuming, labor-intensive
  - Expensive test media and reagents

# Introduction

1996: identification of Gram (-) and Gram (+) bacteria taken directly from culture using MALDI-TOF MS.

Nat Biotechnol. 1996 Nov;14(11):1584-8.

## **The rapid identification of intact microorganisms using mass spectrometry.**

Claydon MA, Davey SN, Edwards-Jones V, Gordon DB.

Department of Biological Sciences, Manchester Metropolitan University, UK.

### **Abstract**

Antibiotic-resistant strains of bacteria continue to emerge, increasing the need for their fast and accurate identification. Matrix-assisted laser desorption ionization time of flight mass spectrometry (MALDI-TOF-MS), has become a prominent technique in biological mass spectrometry. We report the application of MALDI-TOF-MS for the identification of intact Gram-negative and Gram-positive microorganisms taken directly from culture. Analysis of bacteria from a single colony is possible, allowing the screening of mixed cultures. Sample preparation is simple and the analysis automated, providing spectra within minutes. The spectra obtained allow identification of microorganisms from different genera, different species, and from different strains of the same species. The procedure provides a unique mass spectral fingerprint of the microorganism, produced from desorbed components of the cell wall. Consistent data were obtained from subcultures grown for 3-day and 6-day periods, from the same cultures 1 day later and from fresh subcultures 2 months later.

# Introduction

2008: first publication using Biotyper strategy (Bruker).

J Clin Microbiol. 2008 Jun;46(6):1946-54. doi: 10.1128/JCM.00157-08. Epub 2008 Apr 9.

## **Evaluation of matrix-assisted laser desorption ionization-time-of-flight mass spectrometry in comparison to 16S rRNA gene sequencing for species identification of nonfermenting bacteria.**

Mellmann A, Cloud J, Maier T, Keckevoet U, Rammingen I, Iwen P, Dunn J, Hall G, Wilson D, Lasala P, Kostrzewa M, Harmsen D.

Institute for Hygiene, University Hospital Muenster, Muenster D-48149, Germany. mellmann@uni-muenster.de


### **Abstract**

Nonfermenting bacteria are ubiquitous environmental opportunists that cause infections in humans, especially compromised patients. Due to their limited biochemical reactivity and different morphotypes, misidentification by classical phenotypic means occurs frequently. Therefore, we evaluated the use of matrix-assisted laser desorption ionization-time-of-flight mass spectrometry (MALDI-TOF MS) for species identification. By using 248 nonfermenting culture collection strains composed of 37 genera most relevant to human infections, a reference database was established for MALDI-TOF MS-based species identification according to the manufacturer's recommendations for microflex measurement and MALDI BioTyper software (Bruker Daltonik GmbH, Leipzig, Germany), i.e., by using a mass range of 2,000 to 20,000 Da and a new pattern-matching algorithm. To evaluate the database, 80 blind-coded clinical nonfermenting bacterial strains were analyzed. As a reference method for species designation, partial 16S rRNA gene sequencing was applied. By 16S rRNA gene sequencing, 57 of the 80 isolates produced a unique species identification ( $\geq 99\%$  sequence similarity); 11 further isolates gave ambiguous results at this threshold and were rated as identified to the genus level only. Ten isolates were identified to the genus level ( $\geq 97\%$  similarity); and two isolates had similarity values below this threshold, were counted as not identified, and were excluded from further analysis. MALDI-TOF MS identified 67 of the 78 isolates (85.9%) included, in agreement with the results of the reference method; 9 were misidentified and 2 were unidentified. The identities of 10 randomly selected strains were 100% correct when three different mass spectrometers and four different cultivation media were used. Thus, MALDI-TOF MS-based species identification of nonfermenting bacteria provided accurate and reproducible results within 10 min without any substantial costs for consumables.

# Introduction

## What is the Biotyper strategy?


- Identifies microorganisms using MALDI-TOF MS to measure a unique molecular fingerprint of an organism.
- Measures highly abundant proteins that are found in all microorganisms.
- Characteristic patterns of these highly abundant proteins are used to reliably and accurately identify a particular microorganism down to the species level by matching the respective pattern with an extensive open database.





# Sample preparation

Culture stored at RT or 4°C



## Method 1a:

direct transfer of one CFU

Dead easy! The matrix kills the bacteria.

For non pathogenic bacteria (the goodies)

## Method 1b:

Method 1 + 1uL FA

In case method 1 fails...

## Method 2:

EtOH/FA extraction

FA lyses cell wall and releases proteins.

For non sporulating bacteria (the badies)

## Method 3:

TFA extraction

TFA, stronger acid than FA, breaks off the spores.

For sporulating bacteria (the nasties!)

**NB:** method 2 is claimed to be the most efficient. If safe lab practices are observed, method 1 can be used for any type of bacteria.

# Sample spotting



Spot CFU or supernatant  
onto target plate.

↓  
Air dry.

↓  
cover with 2  $\mu$ L HCCA matrix solution.

↓  
Air dry.

↓  
MALDI-TOF analysis




HCCA:  $\alpha$ -cyano-4-hydroxycinnamic acid

MTP 384 polished steel target plate


# MALDI-TOF analysis

- UltrafleXtreme MALDI-TOF/TOF MS (Bruker)
- 10  $\mu\text{m}$  laser diameter for higher resolution
  - 1 kHz speed





# Biotyper workflow



1/ Create a new classification project

# Biotyper workflow

**MALDI Biotyper Realtime Classification Wizard**

**Definition of Project**

Please create a new project for your classification run. The project will record all results generated here. (To continue an existing project, please select the appropriate project name.)

Project Name: 2013-11-08\_BTS Open Tree...

Description: bacterial test standard, 5 reps, Delphine.

Creation Date: 11/8/2013 1:35:07 PM

Creator: Administrator@MALDI-PC

New... Delete Import... View Results...

Next >> Cancel Help

**New Project**

Name: 2013-11-08\_BTS

Creator: Administrator@MALDI-PC

Description: bacterial test standard, 5 reps, Delphine.


OK Cancel Help

# Biotyper workflow

**MALDI Biotyper Realtime Classification Wizard**

**Analyte Placement**

Please specify the target positions for your analytes by drawing a rectangle, clicking on row/column names or directly on the appropriate spots. Analytes are inserted using the Insert key or with Add Analytes from the context menu. Available target spots are shown in hollow. Spots containing analytes of the current project are white (not yet measured) or green (already measured). Please fill in the ID column if it is empty.



2/ List the spotted analytes

|   | Position | Chip | Creation Date | Name [optional] | ID | Description [optional] |
|---|----------|------|---------------|-----------------|----|------------------------|
|   | A1       | 0    | -             |                 |    |                        |
|   | A2       | 0    | -             |                 |    |                        |
|   | A3       | 0    | -             |                 |    |                        |
|   | A4       | 0    | -             |                 |    |                        |
| ▶ | A5       | 0    | -             |                 |    |                        |


Validation Position:

# Biotyper workflow

**MALDI Biotyper Realtime Classification Wizard**

### Analyte Placement

Please specify the target positions for your analytes by drawing a rectangle, clicking on row/column names or directly on the appropriate spots. Analytes are inserted using the Insert key or with Add Analytes from the context menu. Available target spots are shown in hollow. Spots containing analytes of the current project are white (not yet measured) or green (already measured). Please fill in the ID column if it is empty.



|  | Position | Chip | Creation Date | Name [optional] | ID    | Description [optional] |
|--|----------|------|---------------|-----------------|-------|------------------------|
|  | A1       | 0    | -             |                 | BTS-1 |                        |
|  | A2       | 0    | -             |                 | BTS-2 |                        |
|  | A3       | 0    | -             |                 | BTS-3 |                        |
|  | A4       | 0    | -             |                 | BTS-4 |                        |
|  | A5       | 0    | -             |                 | BTS-5 |                        |

Validation Position:

# Biotyper workflow

**MALDI Biotyper Realtime Classification Wizard**

**Selection of MALDI Biotyper Methods**

Please select the MALDI Biotyper methods for data preprocessing (peak picking) and MSP identification (classification).  
[The methods may be changed using the two buttons on the right side.] Furthermore the MSP Source is defined here,  
specifying which reference patterns are applied in classification.

Preprocessing Method: BioTyper Preprocessing Standard Method ...

MSP Identification Method: MALDI Biotyper MSP Identification Standard Method ...

**MSP Source**

☐ MSPs from Libraries

Listeria  
IVD  
BDAL

☒ MSPs from Taxonomy Trees

Taxonomy  
Bruker Taxonomy  
Projects

Cancel << Back Next >> Finish Help

3/ Select the parameters  
(processing and identification  
methods set by Bruker)



# Biotyper workflow




## DB 4613: complete genus/species list

**Abiotrophia** defectiva  
**Acetobacter** aceti  
Acetobacter pasteurianus  
**Acholeplasma** laidlawii  
**Achromobacter** denitrificans  
Achromobacter insolitus  
Achromobacter piechaudii  
Achromobacter ruhlandii  
Achromobacter sp  
Achromobacter spanius  
Achromobacter xylooxidans  
**Acidaminococcus** fermentans  
Acidaminococcus intestini  
**Acidiphilium** acidophilum  
**Acidovorax** avenae ssp avenae  
Acidovorax defluvi  
Acidovorax delafieldii  
Acidovorax facilis  
Acidovorax konjaci  
Acidovorax temperans  
**Acinetobacter** baumannii  
Acinetobacter baylyi  
Acinetobacter bouvetii  
Acinetobacter calcoaceticus  
Acinetobacter gernerii  
Acinetobacter guillouiae  
Acinetobacter haemolyticus  
Acinetobacter johnsonii  
Acinetobacter junii  
Acinetobacter lwoffii  
Acinetobacter nosocomialis  
Acinetobacter parvus  
Acinetobacter pittii

Brucker Taxonomy database contains  
4613 species of bacteria

# Biotyper workflow



**MALDI Biotyper Realtime Classification Wizard**

**Project Summary**


Please check your project settings. Clicking 'Finish' will start the automatic measurement and classification.

|                                 |   |
|---------------------------------|---|
| Project Name:                   | 2013-11-08_BTS                                    |
| Description:                    | bacterial test standard, 5 reps, Delphine.        |
| Number of Positions to measure: | 5   |
| Validation Position:            | (none)  |
| Preprocessing Method:           | BioTyper Preprocessing Standard Method            |
| MSP Identification Method:      | MALDI Biotyper MSP Identification Standard Method |
| MSP Source:                     | Bruker Taxonomy(Taxonomy Tree)                    |

4/ Finalise and launch automated analyses by clicking on “Finish”

Cancel << Back **Finish** Help

# Biotyper workflow



# Biotyper workflow

**MALDI Biotyper Realtime Classification** Project: 2013-11-08\_BTS

File View Tools Help

1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 21 22 23 24

A B C D E F G H I J K L M N O P

☐ Not occupied  
☐ Prepared  
☒ Aborted  
☒ Measured  
☒ Zeroline spectrum  
☒ Measured, classified green  
☒ Measured, classified yellow  
☒ Measured, classified red  
☒ Zeroline spectrum, not classified

**MSP Info**

|                    |  |
|--------------------|--|
| Comment            | Submitted as Type Strain; Matching hint: |
| Conserve           | False                                    |
| Determined by      | DSM                                      |
| Provided by        | HAM                                      |
| Organism           | Erwinia papayae                          |
| Strain             | DSM 16540T                               |
| Matrix             | HCCA                                     |
| Extraction method  | EtOH-FA                                  |
| Growing conditions | ?  |

**Comment**  
Submitted as Type Strain;  
Matching hint: Species mallotivora / papayae of the genus Erwinia have...

Hide Identified Validation: set verified

| ID    | Name | Position | Chip | Detected Species | Score | Comment  |
|-------|------|----------|------|------------------|-------|--|
| BTS-1 |      | A1       | 0    | no peaks found   |       |  |
| BTS-2 |      | A2       | 0    | Escherichia coli | 2.155 | closely related to Shigella and not definitely distinguishable at the moment |
| BTS-3 |      | A3       | 0    | no peaks found   |       |  |
| BTS-4 |      | A4       | 0    | Escherichia coli | 2.270 | closely related to Shigella and not definitely distinguishable at the moment |
| BTS-5 |      | A5       | 0    | no peaks found   |       |  |

Ready. Server: localhost Port: 8080

# Biotyper workflow

**MALDI Biotyper Realtime Classification** Project: 2013-11-08\_BTS

File View Tools Help

1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 21 22 23 24

A B C D E F G H I J K L M N O P

☐ Not occupied  
☐ Prepared  
☒ Aborted  
☒ Measured  
☒ Zeroline spectrum  
☒ Measured, classified green  
☒ Measured, classified yellow  
☒ Measured, classified red  
☒ Zeroline spectrum, not classified

**MSP Info**

|                    |  |
|--------------------|--|
| Comment            | Matching hint: closely related to Shigella   |
| Conserve           | False  |
| Determined by      | ATCC   |
| Provided by        | THL  |
| Organism           | Escherichia coli   |
| Strain             | ATCC 25922   |
| Matrix             | HCCA   |
| Extraction method  | EtOH-FA  |
| Growing conditions | according to <a href="http://www.lgcpromotech.com">http://www.lgcpromotech.com</a> |

**Comment**  
Matching hint: closely related to Shigella and not definitely distinguishable at the moment

Hide Identified Validation: set verified

| ID    | Name | Position | Chip | Detected Species | Score | Comment  |
|-------|------|----------|------|------------------|-------|--|
| BTS-1 |      | A1       | 0    | no peaks found   |       |  |
| BTS-2 |      | A2       | 0    | Escherichia coli | 2.155 | closely related to Shigella and not definitely distinguishable at the moment |

| Score | Detected Species                       | Comment  | Link                |
|-------|--|--|---------------------|
| 2.155 | Escherichia coli ATCC 25922 THL        | closely related to Shigella and not definitely distinguishable at the moment | <a href="#">562</a> |
| 2.152 | Escherichia coli W3350 MMG             | closely related to Shigella and not definitely distinguishable at the moment | <a href="#">562</a> |
| 2.149 | Escherichia coli MB11464_1 CHB         | closely related to Shigella and not definitely distinguishable at the moment | <a href="#">562</a> |
| 2.090 | Escherichia coli ATCC 35218 CHB        | closely related to Shigella and not definitely distinguishable at the moment | <a href="#">562</a> |
| 2.055 | Escherichia coli DSM 30083T HAM        | closely related to Shigella and not definitely distinguishable at the moment | <a href="#">562</a> |
| 1.995 | Escherichia coli ESBLEA_RSS_1528T CHB  | closely related to Shigella and not definitely distinguishable at the moment | <a href="#">562</a> |
| 1.988 | Escherichia coli ATCC 25922 CHB        | closely related to Shigella and not definitely distinguishable at the moment | <a href="#">562</a> |
| 1.917 | Escherichia fergusonii DSM 13698T HAM  | closely related to Shigella and not definitely distinguishable at the moment | <a href="#">564</a> |
| 1.796 | Escherichia coli RV412_A1_2010_06a LBK | closely related to Shigella and not definitely distinguishable at the moment | <a href="#">562</a> |
| 1.761 | Escherichia coli Nissl VML             | closely related to Shigella and not definitely distinguishable at the moment | <a href="#">562</a> |
| 0.000 | no peaks found                         |  | x                   |

BTS-3

Ready.

Server: localhost Port: 8080



# Biotyper workflow



## Bruker Daltonik MALDI Biotyper Classification Results

### Project Info

|                             |  |
|-----------------------------|--|
| Project Name:               | 2013-11-08_BTS                             |
| Project Description:        | bacterial test standard, 5 reps, Delphine. |
| Project Owner:              | Administrator@MALDI-PC                     |
| Project Creation Date/Time: | 2013-11-08T13:35:07.156                    |
| Project Analyte Count:      | 5  |
| Project Type:               | Development                                |
| Validation:                 | not present                                |
| Validation Position:        |  |

### Result Overview

| Analyte Name                   | Analyte ID | Organism (best match)            | Score Value | Organism (second best match)     | Score Value |
|--------------------------------|------------|----------------------------------|-------------|----------------------------------|-------------|
| <a href="#">A1</a><br>(-) (C)  | BTS-1      | no peaks found                   | ≤ 0         | no peaks found                   | ≤ 0         |
| <a href="#">A2</a><br>(++) (C) | BTS-2      | <a href="#">Escherichia coli</a> | 2.155       | <a href="#">Escherichia coli</a> | 2.152       |
| <a href="#">A3</a><br>(-) (C)  | BTS-3      | no peaks found                   | ≤ 0         | no peaks found                   | ≤ 0         |
| <a href="#">A4</a><br>(++) (C) | BTS-4      | <a href="#">Escherichia coli</a> | 2.227       | <a href="#">Escherichia coli</a> | 2.227       |
| <a href="#">A5</a><br>(-) (C)  | BTS-5      | no peaks found                   | ≤ 0         | no peaks found                   | ≤ 0         |

# Biotyper workflow




## Meaning of Score Values

| Range           | Description  | Symbols |
|-----------------|--|---------|
| 2.300 ... 3.000 | highly probable species identification                       | (+++)   |
| 2.000 ... 2.299 | secure genus identification, probable species identification | (++)    |
| 1.700 ... 1.999 | probable genus identification                                | (+)     |
| 0.000 ... 1.699 | not reliable identification                                  | (-)     |


## Meaning of Consistency Categories (A - C)

| Category | Description   |
|----------|---|
| A        | <b>Species Consistency:</b> The best match was classified as 'green' (see above). Further 'green' matches are of the same species as the first one. Further 'yellow' matches are at least of the same genus as the first one.           |
| B        | <b>Genus Consistency:</b> The best match was classified as 'green' or 'yellow' (see above). Further 'green' or 'yellow' matches have at least the same genus as the first one. The conditions of species consistency are not fulfilled. |
| C        | <b>No Consistency:</b> Neither species nor genus consistency (Please check for synonyms of names or microbial mixture).   |


# Biotyper workflow



# Biotyper workflow




# Results



16/18 species from the reference set were successfully identified, thus confirming the biochemical assay results.




# Results



Methods 1a (direct transfer), 1b (DT + FA) and 2 (EtOH/FA) are comparable. Method 3 (TFA) gave no identification results.

# Results



Scores are higher for *S. aureus* than for *E. coli*.  
Scores are slightly better when plate are stored at 4C.  
Scores seem to decrease as culture age.

# Conclusions

- Could Biotyper become a routine technique to identify bacteria from animal tissues?
  - MALDI Biotyper reliably identifies bacteria at the species level (16/18 species from the reference set).
  - Method 1 and 2 gives best scores. Method 1 is much simpler and time/cost-efficient.
  - Storage conditions (4°C or RT) and age of culture plate have a negligible effect on results.
  - It's much quicker, easier, cheaper than the traditional methods

YES, it could!

# Future directions

- What's next?
  - Complete Biotyper assessment using the rest of the reference set.
  - Further validate the strategy on “real” samples (from animal –Simone Warner-, or plant –Brendan Rodoni-tissues).
  - Beside bacteria, attempt identification of filamentous fungi and mycobacteria (need to purchase the database).
  - An honour's student might be hired for this.