

BOTTOM-UP AND TOP-DOWN ANALYSES OF MILK PROTEINS USING LC-MS/MS

Delphine Vincent¹, Aaron Elkins¹, Mark Condina², Benjamin G Cocks^{1,3}, Simone Rochfort^{1,3}


1. DEDJTR, Biosciences Research Division, AgriBio, Bundoora, VIC, Australia
2. Bruker Pty Ltd, Preston, VIC, Australia
3. School of Applied Systems Biology, La Trobe University, Bundoora, VIC, Australia

TOP-DOWN (TD)

TD INTRODUCTION


The development of robust analytical methods for the separation and quantification of variants of milk proteins is of interest. This work aims at establishing a LC-MS-based method to study intact major proteins (caseins and whey) from cow's skim milk. High resolution spectra subsequently allowed accurate mass determination and top-down sequencing (TDS) of major proteins. This LC-MS workflow enables a high-throughput profiling of abundant proteins of cow's milk with minimal sample preparation.

Figure TD1: Experimental workflow.




Skim milk samples were prepared according to Bohe *et al.* (J. Agric. Food Chem. 1998, 46, 458-463).

Figure TD2: Ionisation efficiency (A-B), and LOD and LOQ of milk protein standards (C).




The LC-MS analysis of milk protein standards showed little co-elution (A, B). The ionisation efficiency varied across standards (A), yet it was not concentration-dependent (B). The limit of detection (LOD) ranged from 0.28 (bLG) to 1.35 (aCN) mg/mL and the limit of quantitation (LOQ) ranged from 0.86 (bLG) to 4.09 (aCN) mg/mL (C).

Figure TD3: Mass spectrum deconvolution of aLA elution peak (A) and list of accurate masses from milk standard (B).




Mass spectra were isotopically resolved (A) which allowed for accurate monoisotopic mass determination by MaxEnt deconvolution followed by SNAP annotation of the standards. Each standard eluted several peaks which were successfully deconvoluted into variants and quantified (B).

Figure TD4: LC separation of milk samples.



Milk samples were reproducible and produce 14 chromatographic peaks (A) containing 28 overlapping proteins (B) whose abundance varied across breeds.

Figure TD5: MS/MS analysis of aS1CN.



aS1CN amino acid sequence (Mature protein without peptide signal). In red are shown sequencing results. RPKHPKIQGLQAEVLNLRFFVFPPEVFGKEKVNLSKIDGSESTEDQMEDIKOMEASISSSEIVPNSVEQKHQEDVPSEYLGYLEQLRLKKYKVPQLEVPNSAEERLHSMKEGHIHQAKPEMIGVNOELAYFYPLRFQYQLDAYPSGAYVYPLQTDYTPDSFSDIPNPGENSEKTTMPLW

MS/MS fragmentation using in-source cell-induced-collision (ISCID) was successfully obtained on standards, such as aCN and milk proteins. Resulting a-, y- and b- ion series obtained from ISCID spectra were annotated against the AA sequence using BioTools software (Bruker), as illustrated here for aCN.

TD CONCLUSIONS

In this top-down study, we were able to perform quantitative MS and MS/MS analyses of proteins from Sigma standards and cow's milk samples. MS spectra were of high resolution and accurate protein masses were obtained, thus identifying milk protein isoforms. Top-down sequencing results achieved partial amino acid coverage of C- and N-termini for protein isoform characterisation. This method will be applied in a high-throughput fashion to hundreds of milk samples. Protein variants and their quantities will be related to milk properties and cow phenotype/genotype, and incorporated to dairy cow breeding programs.

ABSTRACT

Cow's milk is a complex fluid whose proteome displays a diverse set of proteins of high abundance such as caseins (CNs) and medium to low abundance whey proteins such as beta-lactoglobulin (bLG), alpha-lactalbumin (aLA), lactoferrin, immunoglobulins, glycoproteins, peptide hormones and enzymes. CNs represent 80% of milk proteins and possess important nutritional and functional properties, whilst bLG (largely unknown function) and aLA (involved in lactose biosynthesis) represents approximately 10 and 4%, respectively. As genetic variants of CNs and whey proteins and post-translational modifications (PTMs) influence many properties of milk that are essential to the dairy industry.


This work had two purposes:

- 1/ Optimising a bottom-up (BU) proteomics procedure to identify as many proteins as possible in a high-throughput fashion.
 - 2/ Establishing a top-down (TD) proteomics method to study intact major proteins, along with the identification of post-translational modifications (PTMs).
- Applying these bottom-up and top-down methods on bulk milk samples from Jersey and Holstein cows, we demonstrate the complementarity of both approaches.

BU INTRODUCTION


A sample preparation method that enables high reproducibility and throughput is key in reliably identifying proteins present in milk. Using skim milk samples from Jersey and Holstein cows, we compared three extraction procedures which have not previously been applied to samples of cows' milk. Method A (urea) involved a simple dilution of the milk in a urea-based buffer, method B (TCA/acetone) involved a TCA/acetone precipitation and method C (methanol/chloroform) involved a tri-phasic partition method in chloroform/methanol solution. Protein assays, SDS-PAGE profiling, and trypsin digestion followed by nanoHPLC-electrospray ionisation-tandem mass spectrometry (nLC-ESI-MS/MS) analyses were performed to assess their efficiency.

Figure BU1: Experimental workflow.




Technical replicates were used during the extraction, digestion, as well as randomised injection steps to assess the reproducibility.

Figure BU2: Comparison of TICs (A), SDS-PAGE patterns (B), protein amount per mL of milk (C), and number of protein accessions per sample (D).




The 3 methods produced distinct TICs (A), with method C displaying peaks with higher resolution than methods A and B. SDS-PAGE patterns (B) were similar from one extraction method to another, albeit extracts C displayed the best resolution with the sharpest bands. In particular extracts C were the only ones consistently resolving the very intense 24-26 kD band corresponding to alpha-caseins, and were therefore more comparable to skim milk profiles than extracts A and B. Method C yielded the highest protein concentrations (C) substantiated by the highest recovery rate, followed by method A, while method B resulted in the lowest concentrations particularly for Holstein breed. Extracts C generated less unique IDs (D) than extracts A, albeit more than B.

Figure BU3: PCA plots along PC1 and PC2 (A), and PC1 and PC7 (B).



Within each method, all replicates clustered together whether it be at the extraction, digestion or injection levels. Within methods, cow breeds did not cluster together.

Figure BU4: Venn diagram of the number of unique protein accessions and Gene Ontology (GO) classification of known proteins per extraction method.



As different as methods A, B and, C were from each other, they recovered the same types of proteins from skim milk samples. In total, 186 different protein accessions were identified across all methods. Identities common to all three sets of extracts include major milk proteins: caseins (aS1, aS2, b, and k forms), lactoferrin, albumin, bLG, aLA, complement C3 and butyrophilin. Yet proteins present in low abundance in milk were also identified, such as enzymes and minor glycoproteins, as well as many immunoglobulins (Igs), antibodies and antigens. All considered, GO classifications were very similar across methods, with method B generally displaying the smallest number of proteins per category. As expected the most prominent protein category in the "Cellular Component" classification was the "extracellular region" as most milk proteins are secreted.

BU CONCLUSIONS

Overall 186 unique accessions, major and minor proteins, were identified across methods. Method C yielded the best resolved SDS-patterns and highest protein recovery rates, method A yielded the greatest number of accessions, and, of the three procedures, method B was the least compatible of all. Our results also highlighted breed differences between the proteins in milk of Jersey and Holstein cows. This bottom-up method will be applied in a high-throughput fashion to hundreds of milk samples in order to better understand the dynamic nature of the milk proteome, with potential applications for milk quality and processing (e.g. cheese, health).

BOTTOM-UP (BU)